



## APPLICATION NOTE

VPE184J046250

### HPLC Analysis of the three Vitamins Thiamine (B<sub>1</sub>), Niacin (B<sub>3</sub>), and Pyridoxine (B<sub>6</sub>) in Beverages and Concentrates

Dr. Felix Kohl, Teamleader Analytics, DOEHLER GmbH, 64295 Darmstadt, Germany

#### Abstract


This application presents a fast and simple method for the determination of three B vitamins in beverages and concentrates. Ion-pair chromatography in reversed-phase mode is used to enable separation of these highly polar compounds from a standard solution and strawberry syrup. Detection is carried out by UV absorption at wavelengths of 260 nm (NIACIN), 290 nm (PYRIDOXINE), and 246 nm (THIAMINE).

#### Keywords


Water-soluble Vitamins  
Reversed-Phase Ion-Pair Chromatography  
Beverages and Food Analysis

#### Compound Information and Chemical Structures


COMPOUND NAME	CLASSIFICATION
NIACIN (B <sub>3</sub> ) / nicotinic acid	water-soluble vitamin / coenzyme precursor NAD and NADP
PYRIDOXINE (B <sub>6</sub> )	water-soluble vitamin / coenzyme amino acid metabolism
THIAMINE (B <sub>1</sub> )	water-soluble vitamin / coenzyme carbohydrate metabolism



NIACIN (B<sub>3</sub>)



PYRIDOXINE (B<sub>6</sub>)



THIAMINE (B<sub>1</sub>)

Figure 1. Chemical structures of the three water-soluble B vitamins NIACIN (B<sub>3</sub>), PYRIDOXINE (B<sub>6</sub>), and THIAMINE (B<sub>1</sub>).



**VDS optilab**  
Chromatographie  
Technik GmbH

#### Contact Information

Phone: +49 (0) 30 55153901

E-mail: [info@vdsoutilab.de](mailto:info@vdsoutilab.de)

Website: [www.vdsoutilab.de](http://www.vdsoutilab.de)

# VDSpher®


Use a better  
column



## Chromatographic Conditions

COLUMN NAME / PARTICLE SIZE	VDSpher PUR 100 C18-H 5 $\mu\text{m}$	ORDER NUMBER
LENGTH $\times$ INNER DIAMETER	250 $\times$ 4.6 mm	VPE184J046250
SEPARATION MODE	reversed-phase ion-pair chromatography	
MOBILE PHASE	8% acetonitrile and 0.5% triethylamine in 10 mM $\text{KH}_2\text{PO}_4$ at pH 2.8 containing 5 mM 1-heptanesulfonic acid	
ELUTION CONDITIONS	isocratic elution	
FLOW RATE	1.0 $\text{mL} \cdot \text{min}^{-1}$	
INJECTION VOLUME	20 $\mu\text{L}$	
TEMPERATURE	ambient	
DETECTION	UV 260 nm <b>NIACIN</b> , 290 nm <b>PYRIDOXINE</b> , 246 nm <b>THIAMINE</b>	
PRESSURE	120 bar	
HPLC SYSTEM	Perkin Elmer Series 200 System Thermo Scientific Accela PDA Detector	
SAMPLE PREPARATION	N/A	

## Chromatographic Analysis



**Figure 1.** LEFT: HPLC chromatogram of a standard solution containing the three water-soluble B vitamins detected at 260 nm (**NIACIN**), 290 nm (**PYRIDOXINE**), and 246 nm (**THIAMINE**). RIGHT: HPLC chromatogram of a strawberry syrup sample containing two of three water-soluble B vitamins detected at 260 nm (**NIACIN**), and 290 nm (**PYRIDOXINE**).

## References

Y. Maeda *et al.*, *J. Assoc. Off. Anal. Chem.* **1989**, 72, 244 – 247.



**VDS optilab**  
Chromatographie  
Technik GmbH

### Contact Information

Phone: +49 (0) 30 55153901

E-mail: [info@vdsoptilab.de](mailto:info@vdsoptilab.de)

Website: [www.vdsoptilab.de](http://www.vdsoptilab.de)