

SEPARA® ONE-STEP-FILTER-VIALS

Easy to use - efficient and cost-saving!

Save time and money in the sample preparation process with the Separa® Filter Vials.

Separa® is a one-step sample preparation device that allows a quick, effective, and simple to use process integrating a self-sampling system, membrane, cap, and lid with pre-cut septa.

The liquid is placed in the outer vial and dispensed by pressing the inner insert so that it is pushed through the membrane. The sample is then filtered and Separa® is ready to be loaded into the autosampler.

Einfach in der Anwendung - effizient und kostensparend!

Sepa® Filter Vials - die Ein-Schritt-Probenvorbereitung, die Zeit und Geld spart!

Die in fünf unterschiedlichen Membranausführungen erhältlichen Separa® Filter-Vials sind mit vorgeschlitzten Septen ausgestattet und sehr einfach in der Handhabung.

Probenflüssigkeit wird in die äußere Probenflasche gegeben und durch Herunterdrücken des Filtereinsatzes durch die Membran gepresst. Die gefilterte Probenflüssigkeit steht nun im inneren Einsatz des Filter-Vials für die Probenaufgabe durch den Autosampler bereit.



Part No. Art.-Nr.	Description Artikelbeschreibung	Amount Menge	Price Preis
FV-32PT020S04-1	32 mm PTFE FILTER Vials 0.2 µm, Pink, 100/pk	100	205.00 €
FV-32PT020S04-5	32 mm PTFE FILTER Vials 0.2 µm, Pink, 500/pk	500	894.00 €
FV-32PT045S05-1	32 mm PTFE FILTER Vials 0.45 µm, Red, 100/pk	100	198.00 €
FV-32PT045S05-5	32 mm PTFE FILTER Vials 0.45 µm, Red, 500/pk	500	869.00 €
FV-32PV020S04-1	32 mm PVDF FILTER Vials 0.2 µm, Yellow, 100/pk	100	205.00 €
FV-32PV020S04-5	32 mm PVDF FILTER Vials 0.2 µm, Yellow, 500/pk	500	894.00 €
FV-32PV045S05-1	32 mm PVDF FILTER Vials 0.45 µm, Orange, 100/pk	100	205.00 €
FV-32PV045S05-5	32 mm PVDF FILTER Vials 0.45 µm, Orange, 500/pk	500	894.00 €
FV-32PS020S04-1	32 mm PES FILTER Vials 0.2 µm, Light Green, 100/pk	100	205.00 €
FV-32PS020S04-5	32 mm PES FILTER Vials 0.2 µm, Light Green, 500/pk	500	894.00 €
FV-32PS045S05-1	32 mm PES FILTER Vials 0.45 µm, Dark Green, 100/pk	100	205.00 €
FV-32PS045S05-5	32 mm PES FILTER Vials 0.45 µm, Dark Green, 500/pk	500	894.00 €
FV-32NL020S04-1	32 mm Nylon FILTER Vials 0.2 µm, Light Blue, 100/pk	100	205.00 €
FV-32NL020S04-5	32 mm Nylon FILTER Vials 0.2 µm, Light Blue, 500/pk	500	894.00 €
FV-32NL045S05-1	32 mm Nylon FILTER Vials 0.45 µm, Blue, 100/pk	100	205.00 €
FV-32NL045S05-5	32 mm Nylon FILTER Vials 0.45 µm, Blue, 500/pk	500	894.00 €
FV-32RC020S04-1	32 mm RC FILTER Vials 0.2 µm, Gray 100/pk	100	205.00 €
FV-32RC020S04-5	32 mm RC FILTER Vials 0.2 µm, Gray 500/pk	500	894.00 €
FV-32RC045S05-1	32 mm RC FILTER Vials 0.45 µm, Black, 100/pk	100	205.00 €
FV-32RC045S05-5	32 mm RC FILTER Vials 0.45 µm, Black, 500/pk	500	894.00 €

01_2023

Bei größerem Bedarf fordern Sie bitte Ihr individuelles Angebot an. Gerne senden wir Ihnen Filter-Vials zum Test.